



ADENNA®

A hospeco BRANDS GROUP Company

Self Sealing Sterilization Pouch

Product Specification Sheet

SKUs	SPH2338, SPH2350, SPH2390, SPH2810, SPH3510, SPH3565, SPH5311, SPH5385, SPH7514, SPH1016, SPH1219	
Film Thickness	52um ±2% (PET+CPP film)	
Paper	Weight	60 g/m ² ±2% (Medical grade paper)
	Air Permeance	≤60 (Sec/100ml)
	Water Absorbency(60"s)	≤30g/m ²
Peel Strength	before Sterilization	≥ 1.5 (N/15mm)
	after Sterilization	≥ 1.5 (N/15mm)
Width of Sealing Edge	≥ 10mm	
Fiber in Peeling	no	
Sealing Method	Self adhesive	
Sterilization Method	Steam: sterilize for 20 minutes at temperature 121°C sterilize for 3 minutes at temperature 134°C EO: temperature - 30-70°C; humidity - 30-80RH; gas thickness: 400-800mg/L	
Storage Condition	Please store at temperature 25°C± 5°C, humidity 50% ± 5% RH, to keep sure the sterilization effect is not affected.	
Opening Method	Please start opening at the thumb notch, use average power for both sides, tear off slowly and softly till 1/2 or 2/3 point of the pouch, then take out the object.	
Date of Validity	It can be kept for 3 years before sterilization; and 6 months after sterilization.	

Note: Please follow the instruction sheet to store and use the material, to make sure the best sterilization effect.